

MEYRA



FLASH

+ **SUPER FLEXIBLE**
A SEAT THAT GROWS WITH YOU



Crash tested for safety
Complies with international
ISO safety standards for
crash safety

- +** **Highly adjustable growth seating**
- +** **Adjustable seat angles**
- +** **Adjustable wheel positions for balance**
- +** **Adjustable wheel camber**
- +** **High quality German design and manufacture**
- +** **Crash tested—user can be transported in their chair by bus and maxi taxi**

Compact & nimble

for small kids to big teens

ADAPTABLE TO EACH CHILD'S NEEDS

Designed for growth, the Meyra Flash is a highly adjustable manual wheelchair with lots of options to customise.

With the Flash, you can have different floor to seat heights, front and rear; different axle positions for changing seat angle/dump; different rim, wheel and tyre combinations; user and attendant brake options; spoke guards; anti-tip wheels—and a variety of frame colours!



Detachable anatomic back shell made of aluminium



Cool design. Low weight. Lightning quick and agile.



Infinitely adjustable lower leg length, continuous from 200 mm



Super flexible—a seat that grows with you

TECHNICAL DATA

FLASH

	midi frame	maxi frame
Seat width	24 - 30 cms	24 - 40 cms
Seat depth	24 - 30 cms	32 - 40 cms
Seat height - front/rear	36-44/39-47cms	39-48/39-47cms
Backrest height	25 - 40 cms	25 - 40 cms
Backrest angle adjustment	75 - 105 degrees	75 - 105 degrees
Legrest length	20 - 29 cms	25 - 32/33 - 44cms
Overall length - with legrests	Approx. 67 - 75cms	Approx. 86-93cms
Wheel size rear	55cms/22" 60cms/24"	
Chair weight	Approx. 10 kgs	
Maximum client weight	75 kgs	
Compliance testing	Complies with ISO7176 strength & durability standards Tested and compliant with ISO7176-19 crash test requirements	

Dimension tolerances $\pm 10 \text{ mm} \pm 2^\circ$, Data subject to change

MEYRA
AUSTRALIA

MEYRA Australia
5 Boulder Rd Malaga
Western Australia 6090

Tel (08) 9248 3560
Fax (08) 9209 1999

sales@meyraaustralia.com.au
www.meyraaustralia.com.au