

MEYRA



MC-S

+ COLOURFUL
AND COMFORTABLE



Crash tested for safety
Complies with international
ISO safety standards for
crash safety

- +** *Strong, compact, safe and reliable*
- +** *Adjustable—seat frame grows with the child, customisable with different seating systems*
- +** *Wide range of frame colours & accessories*
- +** *Easy to drive and control with standard and optional attendant controls*
- +** *High quality German design and manufacture*
- +** *Crash tested—can be transported by bus and maxi taxi*

A little gem

COLOURFUL, ADJUSTABLE AND SAFE

A strong, safe, durable and reliable power wheelchair is essential for kids.

The Meyra MC-S power wheelchair is all this and more! It has a growth-accommodating seating system, it's easy to drive, has excellent manoeuvrability, easy access; it's also easy to clean, maintain and manage.

Make the MC-S part of your family!



Central leg support, 90° knee angle



200 mm seat lift



Tilt-in-space



5 frame colours

TECHNICAL DATA

MCS

Seat width	30 - 44 cms
Seat depth	30 - 40 cms
Backrest height	39 - 44 cms
Legrest length	25 - 35 cms
Armrest height	18 - 25 cms
Seat height - front	42 - 47 cms
Overall length - with legrests	Approx. 97 cms
Overall length - without legrests	Approx. 79 cms
Overall width	58 cms
Wheel size front	190 x 50mm/225 x 70mm
Wheel size rear	320 X 60MM
Maximum obstacle climbing height	Approx. 60mm

*distance will be affected by client weight, terrain and battery condition
Tolerance ± 10 mm
Data subject to change

Turning radius	Approx. 84 cms
Driving range	Approx. 35 to 40 kms with 73 ah batteries*
Chair weight	Approx. 89 kgs, no batteries, armrests, legrests
Maximum client weight	75 kgs
Motor power rating	300 watt German motors
Maximum gradient	15 degrees
Mechanical seat tilt	0 to 10 degrees
Power seat tilt	0 to 22 degrees, tilt in space
Power seat raise & tilt	20 cms raise, 30 degree tilt
Mechanical backrest	-10 to + 30 degrees
Power backrest recline	-5 to + 35 degrees
Compliance testing	Complies with ISO7176 strength & durability standards Tested and compliant with ISO7176-19 crash test requirements

MEYRA
AUSTRALIA

MEYRA Australia
5 Boulder Rd Malaga
Western Australia 6090

Tel (08) 9248 3560
Fax (08) 9209 1999

sales@meyraaustralia.com.au
www.meyraaustralia.com.au