

MEYRA

ORBIT



NIMBLE
COLOURFUL
RELIABLE



Crash tested



Precise manoeuvring

Manoeuvre around obstacles or turn on the spot



Ideal for confined spaces

Low seat height and overall width



Colours how you want them

Pick your colours—drive wheel rims, frame accents, clothing guard



Crash tested for safety

Complies with international ISO safety standards for crash safety


OUTSTANDINGLY NIMBLE

THE AGILE & COLOURFUL MID-WHEELER

Meyra's Orbit will draw you in with its agile mid-wheel drive configuration and striking colour accent choices. It's no lightweight either, with a maximum user weight of 160 kgs. Built with Meyra's world-famous German engineering prowess, this is the mid-wheel drive you've been waiting for.



Manoeuvre around obstacles or turn on the spot



New ErgoSeat cover with grey seams



Colours - wheel rim, frame accent & guard accent



Single or double support castor

TECHNICAL DATA

Seat width ¹	380 - 650
Seat depth ⁴	400 - 560 / 400 - 560 / 480 - 530
Seat height front ⁶	430, 465, 500 / 440, 475, 510
Back height ⁴	450 - 500 / 530 - 580 / 640
Lower leg length ²	280 - 450
Armrest height ²	240 - 350
Length with footrests approx. ⁷	1,180 - 1,250 / 1,060 - 1,150
Length without footrests	910
Width, ready to move min. ⁸	600/600/630/680 with 12" 640/640/640/680 with 14"
Height without headrest ³	980 - 1,060
Transport height min.	630
Obstacle height	60
Turning radius (with central legrest)	700
Turning radius (with single legrest)	800
Ground clearance	70

Min. transport weight empty	120 - 165
Max. user weight incl. additional load	120 / 160 ⁵
Speed in km/h	6 / 10
Motor power W	2 x 180 / 2x 350
Range with battery 50Ah (20h) max. in km	25 - 30
Range with battery 73Ah (20h) max. in km	35 - 40
Permitted uphill / downhill gradient	15%
Steering type	Indirect
Seat tilt, electrical	0° to 30° / 0° to 50°
Seat tilt, mechanical	1°, 4°, 7°, 10°
Stand-up aid (0° to -20°) with electrical tilt, lift 200 mm	0° to 45°
Back angle, electrical	10° to 50°
Lift	300

1) see category Seat width 01
2) to seat pan without cushion
3) with adjustable back
4) seating systems:
Standard / ErgoSeat / Recaro

5) with 350W motors
6) measured at 0° seat angle, 12"/14"
7) Single/central legrests
8) Code 38/43/48/53, with 350W
+10 mm total width

Dimensions in mm, weights in kg, wheels in inches, speed in km/h
Dimension tolerances ± 15 mm, +/-2°, data subject to changes.

MEYRA
AUSTRALIA

MEYRA Australia
5 Boulder Rd Malaga
Western Australia 6090



Tel (08) 9248 3560
Fax (08) 9209 1999



sales@meyraaustralia.com.au
www.meyraaustralia.com.au